

ABX **Micros ES** 60  
ABX **Micros ESV** 60  
**Micros Care ST**  
**sci1 Vet abc Plus<sup>+</sup>**

## Printer list

Ref.: RAL095E

## **Drivers available by default on the instrument**

- deskjet
- djet500
- laserjet
- lj5gray
- lj5mono
- ljet2p
- ljet3
- ljet3d
- ljet4
- ljet4d
- ljetplus

## **New drivers provided in the update**

- bbox
- bit
- bitcmyk
- bitrgb
- bj10e
- bj200
- bjc600
- bjc800
- cdeskjet
- cdj550
- cdjcolor
- cdjmono
- cljet5
- cljet5c
- devicen
- epson
- epsonc
- eps9mid
- eps9high
- epswrite
- faxg3
- faxg32d
- faxg4
- ibmpro
- lx300
- ml-320
- okiibm
- pj
- pjl
- pjl300
- psdcmyk
- psdrgb
- psgray
- psmono
- psrgb
- pswrite
- pxcolor
- pxlmono
- spotcmyk
- xcf

## Printers validated by HORIBA Medical

The compatibility of the printers listed below has been tested. HORIBA Medical can guarantee the functioning of ABX Micros ES60 OT/CT, ABX Micros ESV60, Micros Care ST and scil Vet abc Plus<sup>+</sup> instruments with these printers:

Printer	HORIBA Medical part number	Driver
Epson EPL 5800 (laser)	Obsolete	ljetplus, ljet4d, ljet4, ljet3d, ljet3, ljet2p, lj5mono, lj5gray, laserjet, djet500, deskjet
Epson EPL 5900 (laser)	Obsolete	ljetplus, ljet4d, ljet4, ljet3d, ljet3, ljet2p, lj5mono, lj5gray, laserjet, djet500, deskjet
Epson EPL 6200 110V (laser) Epson EPL 6200 220V (laser)	XAA635A XAA541B	ljet4
Epson LX 300+ 110V (USB) Epson LX 300+ 220V (USB)	XAA476A XAA473A	lx-300
OKI B430d 110V OKI B430d 220V	XAA627A XAA626A	ljet4
OKI B4200 (laser)	Obsolete	ljetplus, ljet4d, ljet4, ljet2p, lj5mono, lj5gray, laserjet, djet500, deskjet
OKI B4350 (laser)	Obsolete	ljetplus, ljet4d, ljet4, ljet3d, ljet3, ljet2p, lj5mono, lj5gray, laserjet, djet500, deskjet
OKI B4600 110V OKI B4600 220V	XAA594A XAA593A	ljet4, pxlmono
OKI Microline 320 turbo	/	okiibm

Other printers\* (not tested by HORIBA Medical) can be used with the drivers integrated in the instruments. They are sorted by printer trademark:

- ["Brother, page 15"](#)
- ["Canon, page 17"](#)
- ["Epson, page 20"](#)
- ["Fujitsu, page 21"](#)
- ["Generic, page 22"](#)
- ["HP, page 23"](#)
- ["IBM, page 26"](#)
- ["Kyocera, page 27"](#)
- ["Lexmark, page 31"](#)
- ["Minolta, page 32"](#)
- ["NEC, page 33"](#)
- ["Oce, page 34"](#)
- ["Okidata, page 35"](#)
- ["Panasonic, page 37"](#)
- ["Samsung, page 38"](#)
- ["Sharp, page 40"](#)
- ["Star, page 43"](#)
- ["Toshiba, page 44"](#)
- ["Xerox, page 45"](#)
- ["Others, page 47"](#)

\* This list has been established in June 2011 according to the information provided in the Linux website ([http://www.linuxprinting.org/driver\\_list.cgi](http://www.linuxprinting.org/driver_list.cgi)).



To know the new printers added, please refer to ["New Printers, page 5"](#). These new printers are also highlighted in blue in printer trademark tables.



**Caution: some printers may require the purchase of USB cord.**

## New Printers

Printer	Driver
Anitech M24	epson
Apollo P-2200	cdj550
Apollo P-2250	cdj550
Brother HL-2030	lj5gray, ljet4
Brother HL-2035	lj5gray, ljet4
Brother HL-2040	ljet4
Brother HL-2140	lj5gray, ljet4, pxlmono
Brother HL-2170W	ljet4, ljet4d
Brother HL-4040CN	pxlcolor
Brother HL-4050CDN	pxlcolor
Brother HL-4070CDW	pxlcolor
Brother MC-3000	epson
Brother MFC-7150C	epsonc
Brother MFC-9100c	epsonc
Brother MFX-9840CDW	pxlcolor
Canon BJ-5	bj10e
Canon BJ-10e	bj10e
Canon BJ-20	bj10e
Canon BJ-30	bj200
Canon BJ-100	bj200
Canon BJ-200	bj200
Canon BJC-210	bjc600
Canon BJC-240	bjc600
Canon BJC-250	bjc600
Canon BJ-330	bj200
Canon BJC-50	bjc600
Canon BJC-55	bjc600
Canon BJC-70	bj200, bjc600
Canon BJC-80	bjc600
Canon BJC-85	bjc600
Canon BJC-210	bj200
Canon BJC-210SP	bj200
Canon BJC-240	bj200
Canon BJC-250	bj200
Canon BJC-600	bjc600
Canon BJC-610	bjc600
Canon BJC-620	bjc600
Canon BJC-800	bjc800
Canon BJC-1000	bj200, bjc600
Canon BJC-2000	bjc600
Canon BJC-2010	bjc600
Canon BJC-2100	bjc600
Canon BJC-2110	bjc600

**Printer list**  
New Printers

Printer	Driver
Canon BJC-3000	bjc600
Canon BJC-4000	bjc600
Canon BJC-4100	bjc600
Canon BJC-4200	bjc600
Canon BJC-4300	bjc600
Canon BJC-4310SP	bjc600
Canon BJC-4400	bjc600
Canon BJC-4550	bj200, bjc800
Canon BJC-7000	bjc800
Canon BJC-7100	bjc800
Canon BJC-4550	bjc800
Canon i450	bjc800
Canon imageRunner 1023	pxlmono
Canon imageRunner 1023iF	pxlmono
Canon imageRunner 1023N	pxlmono
Canon imageRunner 2016	pxlmono
Canon imageRunner 2016i	pxlmono
Canon imageRunner 2018	pxlmono
Canon imageRunner 2020	pxlmono
Canon imageRunner 2020i	pxlmono
Canon imageRunner 2022	pxlmono
Canon imageRunner 2230	pxlmono
Canon imageRunner 2270	pxlmono
Canon imageRunner 3025	pxlmono
Canon imageRunner 3225	pxlmono
Canon imageRunner 7086	pxlmono
Canon imageRunner 7095	pxlmono
Canon imageRunner 7105	pxlmono
Canon imageRunner 8070	pxlmono
Canon imageRunner C2550	pxlcolor
Canon imageRunner C3200	cljet5
Canon imageRunner C5058	pxlcolor
Canon imageRunner C5068	pxlcolor
Canon imageRunner C5180	cljet5
Canon imageRunner C5185	pxlcolor
Canon imageRunner C5185i	pxlcolor
Canon LBP-3460	pxlmono
Canon LBP-5360	cljet5
Canon LBP-5960	cljet5
Canon LBP-5970	pxlmono
Canon LBP-5975	pxlmono
Canon MultiPASS C2500	bjc600
Canon MultiPASS C3000	bjc600
Canon MultiPASS C3500	bjc600
Canon MultiPASS C5000	bjc600

**Printer list**  
New Printers

Printer	Driver
Canon MultiPASS C5500	bjc600
Canon S100	bjc600
Canon S400	bjc600
Dell A920	deskjet
Dell 3100cn	pxlcolor
Deskjet 3050 all-in-one J610a series	deskjet
Epson AcuLaser C1900 PS3	ljet4
Epson AcuLaser C2000 PS3	lj5gray, ljet4
Epson AcuLaser C4000 PS3	ljet4
Epson AcuLaser C8600 PS3	lj5gray, ljet4, jlet4d
Epson Dot Matrix	epson, epsonc, eps9mid, eps9high
Epson EPL-5900 PS3	lj5gray, ljet4, pxlmono
Epson EPL-6100 PS3	lj5gray, ljet4, pxlmono
Epson EPL-N2500 PS3	lj5gray, ljet4, pxlmono
Epson L1000	epson
Epson LQ24	epsonc
Epson LQ-500	epson
Epson LQ-570+	epson
Epson LX-1050	epson
Generic ESC/P Dot Matrix Printer	epson, epsonc, eps9mid, eps9high
Generic IMB-Compatible Dot Matrix Printer	ibmpro
HP Business Inkjet 3000	pxlcolor
HP Color LaserJet 5	cljet5
HP Color LaserJet 5M	cljet5
HP Color LaserJet 5500	pxlcolor
HP Deskjet 310	cdj550
HP Deskjet f340	deskjet
HP Deskjet 560C	cdj550
HP Deskjet 610CL	cdj550
HP Deskjet 612C	cdj550
HP Deskjet 648C	cdj550
HP Deskjet 692C	cdj550
HP Deskjet 694C	cdj550
HP Deskjet 812C	cdj550
HP Deskjet 840C	cdj550
HP Deskjet 842C	cdj550
HP Deskjet 930C	cdj550
HP Deskjet 932C	cdj550
HP Deskjet 950C	cdj550
HP Deskjet 970C	cdj550
HP Deskjet 990C	cdj550
HP Deskjet 1200C	pjl300
HP LaserJet 1320	pxlmono
HP LaserJet 2500	pxlcolor
HP LaserJet 3390	lj5gray, ljet4, ljet4d, pxlmono

**Printer list**  
New Printers

Printer	Driver
HP LaserJet 3392	lj5gray, ljet4, ljet4d, pxlmono
HP LaserJet 4200	pxlmono
HP LaserJet 4300	pxlmono
HP LaserJet 4500	pxlcolor
HP LaserJet 4550	pxlcolor
HP LaserJet 4600	pxlcolor
HP LaserJet 5000	pxlcolor
HP LaserJet 5500	pxlcolor
HP LaserJet 5100	pxlmono
HP LaserJet M1522nf MFP	pxlmono
HP LaserJet M1522n MFP	pxlmono
HP OfficeJet 500	cdj550
HP OfficeJet 600	cdj550
HP OfficeJet G55	cdj550
HP OfficeJet G85	cdj550
HP OfficeJet G95	cdj550
HP OfficeJet Pro 1150C	cdj550
HP OfficeJet Pro 1170C	cdj550
HP OfficeJet R45	cdj550
HP OfficeJet T45	cdj550
HP OfficeJet T65	cdj550
HP PaintJet	pj
HP PaintJet XL	pjxl
HP PaintJet XL300	pjxl300
HP PhotoSmart P100	pjxl300
HP PhotoSmart P1000	cdj550
HP PhotoSmart P1100	cdj550
IBM Execjet 4072	bj200
IBM ProPrinterII	ibmpro
IBM 4312	laserjet
Infotec 4651	pxlmono
NEC PinWriter P6	epsonc
NEC PinWriter P6 plus	epsonc
Kyocera Ci-1100	cljet5
Kyocera CS-250ci KPD	cljet5
Kyocera CS-300ci KPD	cljet5
Kyocera CS-400ci KPD	cljet5
Kyocera CS-500ci KPD	cljet5
Kyocera CS-C2525E KPD	cljet5
Kyocera CS-C3225E KPD	cljet5
Kyocera CS-C3232E KPD	cljet5
Kyocera CS-C4035E KPD	cljet5
Kyocera FS-C5015N	cljet5
Kyocera FS-C5016N	cljet5
Kyocera FS-C5020N	cljet5

**Printer list**  
New Printers

Printer	Driver
Kyocera FS-C5025N	cljet5
Kyocera FS-C5030N	cljet5
Kyocera FS-C5100DN KPD	cljet5
Kyocera FS-C5200DN KPD	cljet5
Kyocera FS-C5300DN KPD	cljet5
Kyocera FS-C5350DN KPD	cljet5
Kyocera FS-C5400DN KPD	cljet5
Kyocera FS-C8008N	cljet5
Kyocera FS-C8100DN	cljet5
Kyocera FS-C8100DN+KPD	cljet5
Kyocera KM-C830	cljet5
Kyocera KM-C850	cljet5
Kyocera KM-C2520	cljet5
Kyocera KM-C2525E KPD	cljet5
Kyocera KM-C2630	cljet5
Kyocera KM-C3225	cljet5
Kyocera KM-C3225E KPD	cljet5
Kyocera KM-C3232	cljet5
Kyocera KM-C3232E KPD	cljet5
Kyocera KM-C4035E KPD	cljet5
Kyocera TASKalfa 250ci KPD	cljet5
Kyocera TASKalfa 300ci KPD	cljet5
Kyocera TASKalfa 400ci KPD	cljet5
Kyocera TASKalfa 500ci KPD	cljet5
Lexmark 3000	cdj550
Lexmark C720n	pxlcolor
Lexmark E120n	ljet4, ljet4d, pxlmono
Lexmark E230	pxlmono
Lexmark E238	pxlmono
Lexmark E450dn	lj5gray, pxlmono
Lexmark X342n	pxlmono
Lexmark X642e	pxlmono
Oce fx2080	pxlmono
Oce sx 1480	pxlmono
OKI B410	ljet4, ljet4d, pxlmono
OKI B2400	pxlmono
OKI B2400n	pxlmono
OKI B4200	pxlmono
OKI B4250	pxlmono
OKI B4350	ljet4, pxlmono
OKI B2400	pxlmono
OKI B2400n	pxlmono
OKI B4200	pxlmono
OKI B4250	pxlmono
OKI B2400	pxlmono

**Printer list**  
New Printers

Printer	Driver
OKI B2400n	pxlmono
OKI B4200	pxlmono
OKI C5300	cljet5
OKI C5400	cljet5
OKI C6000	pxlcolor
OKI C6100	pxlcolor
OKI C7300	cljet5
OKI C7350	pxlcolor
OKI C7550	pxlcolor
OKI C8800	pxlcolor
OKI C9200	cljet5
OKI C9400	cljet5
OKI C9650	pxlcolor
OKI C9800	pxlcolor
OKI Microline 192+	ibmpro
OKI Microline IBM compatible 9 pin	okiibm
OKI ML 320	okiibm
OKI ML 321	okiibm
OKI ML 380	epson
OKI Okijet 2500	cdj550
Panasonic DP-180	pxlcolor
Panasonic DP-2330	pxlmono
Panasonic DP-4510	pxlmono
Panasonic DP-4520	pxlmono
Panasonic DP-8016P	pxlmono
Panasonic DP-8020E	pxlmono
Panasonic DP-8045	pxlmono
Panasonic DP-8060	pxlmono
Panasonic DP-C213	pxlmono
Panasonic DP-C262	pxlmono
Panasonic DP-C264	pxlmono
Panasonic DP-C265	pxlcolor
Panasonic DP-C305	pxlcolor
Panasonic DP-C322	pxlcolor
Panasonic DP-C354	pxlmono
Panasonic DP-C405	pxlcolor
Panasonic KX-MC6020	pxlcolor
Panasonic KX-P1123	epson
Panasonic KX-P1124	epson
Panasonic KX-P1150	eps9high
Panasonic KX-P1624	epson
Panasonic KX-P2023	epson
Panasonic KX-P2123	epson
Panasonic KX-P2135	epsonc
Panasonic KX-P2150	ibmpro

**Printer list**  
New Printers

Printer	Driver
QMS LPK-100	eps9mid
Samsung CLP-325w	pxlcolor
Samsung CLP-350	pxlcolor
Samsung CLP-550	pxlcolor
Samsung CLP-550N	pxlcolor
Samsung CLP-650	pxlcolor
Samsung CLP-650N	pxlcolor
Samsung CLP-660N	pxlcolor
Samsung CLP-660ND	pxlcolor
Samsung ML-2570	pxlmono
Samsung ML-2571N	pxlmono
Samsung ML-2850D	pxlmono
Samsung ML-2851ND	pxlmono
Samsung ML-3051N	pxlmono
Samsung ML-3051ND	pxlmono
Samsung ML-3470N	pxlmono
Samsung ML-3471ND	pxlmono
Samsung ML-4050N	pxlmono
Samsung ML-4550	pxlmono
Samsung ML-4551N	pxlmono
Samsung ML-4551ND	pxlmono
Samsung ML-6060	pxlmono
Samsung ML-6060N	pxlmono
Samsung ML-6060S	pxlmono
Sharp AR-155FG	pxlmono
Sharp AR-160M	pxlmono
Sharp AR-163FG	pxlmono
Sharp AR-163G	pxlmono
Sharp AR-168D	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-168S	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-200M	pxlmono
Sharp AR-205FG	pxlmono
Sharp AR-205G	pxlmono
Sharp AR-266FP	pxlmono
Sharp AR-266S	pxlmono
Sharp AR-311FP	pxlmono
Sharp AR-311N	pxlmono
Sharp AR-311S	pxlmono
Sharp AR-351FP	pxlmono
Sharp AR-351N	pxlmono
Sharp AR-351S	pxlmono
Sharp AR-451FP	pxlmono
Sharp AR-451N	pxlmono
Sharp AR-451S	pxlmono
Sharp AR-555M	pxlmono

**Printer list**  
New Printers

Printer	Driver
Sharp AR-555S	pxlmono
Sharp AR-625M	pxlmono
Sharp AR-625S	pxlmono
Sharp AR-705M	pxlmono
Sharp AR-705S	pxlmono
Sharp AR-5220	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-BC260	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-BC320	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C170FP	pxlcolor
Sharp AR-C170M	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C172FP	pxlcolor
Sharp AR-C172M	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C260	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C260F	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C260FP	pxlcolor
Sharp AR-C260M	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C260P	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C260S	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C261F	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C261FP	pxlcolor
Sharp AR-C261M	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C261S	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C262FP	pxlcolor
Sharp AR-C262M	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C262S	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-M150	pxlmono
Sharp AR-M155	pxlmono
Sharp AR-M160	pxlmono
Sharp AR-M161	pxlmono
Sharp AR-M162	pxlmono
Sharp AR-M165	pxlmono
Sharp AR-M205	pxlmono
Sharp AR-M206	pxlmono
Sharp AR-M207	pxlmono
Sharp AR-M236	pxlmono
Sharp AR-M237	pxlmono
Sharp AR-M276	pxlmono
Sharp AR-M277	pxlmono
Sharp AR-M351N	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M351U	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M355N	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M355U	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M451N	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M451U	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M455N	lj5gray, ljet4, ljet4d, pxlmono

**Printer list**  
New Printers

Printer	Driver
Sharp AR-M455U	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M550N	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M550U	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M620N	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M620U	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M700N	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M700U	lj5gray, ljet4, ljet4d, pxlmono
Sharp ARP350	lj5gray, ljet4, ljet4d, pxlmono
Sharp ARP450	lj5gray, ljet4, ljet4d, pxlmono
Sharp MX-2300FG	pxlcolor
Sharp MX-2300G	pxlcolor
Sharp MX-2300N	pxlcolor
Sharp MX-2700FG	pxlcolor
Sharp MX-2700G	pxlcolor
Sharp MX-2700N	pxlcolor
Sharp MX-3500FN	pxlcolor
Sharp MX-3500N	pxlcolor
Sharp MX-3501FN	pxlcolor
Sharp MX-3501N	pxlcolor
Sharp MX-4500FN	pxlcolor
Sharp MX-4500N	pxlcolor
Sharp MX-4501FN	pxlcolor
Sharp MX-4501N	pxlcolor
Star LC24-200	epson, epsonc
Star LC 90	epson
Star NL-10	epson
Star NX-1001	eps9mid, eps9high
Toshiba e-Studio 162	pxlmono
Toshiba e-Studio 162d	pxlmono
Toshiba e-Studio 167	pxlmono
Toshiba e-Studio 203l	pxlmono
Toshiba e-Studio 220cp	pxlcolor
Toshiba e-Studio 233	pxlmono
Toshiba e-Studio 237	pxlmono
Toshiba e-Studio 283	pxlmono
Toshiba e-Studio 353	pxlmono
Toshiba e-Studio 450s	pxlmono
Toshiba e-Studio 453	pxlmono
Toshiba e-Studio 500p	pxlmono
Toshiba e-Studio 500s	pxlmono
Toshiba e-Studio 523	pxlmono
Toshiba e-Studio 523t	pxlmono
Toshiba e-Studio 603	pxlmono
Toshiba e-Studio 603t	pxlmono
Toshiba e-Studio 723	pxlmono

**Printer list**  
New Printers

Printer	Driver
Toshiba e-Studio 723t	pxlmono
Toshiba e-Studio 853	pxlmono
Toshiba e-Studio 901	pxlmono
Toshiba e-Studio 1101	pxlmono
Toshiba e-Studio 1351	pxlmono
Toshiba e-Studio 2830c	pxlcolor
Toshiba e-Studio 3530c	pxlcolor
Toshiba e-Studio 4500c	pxlcolor
Toshiba e-Studio 4520c	pxlcolor
Toshiba e-Studio 5500c	pxlcolor
Xerox Document Homecentre	cdj550
Xerox DocuPrint C6	cdj550
Xerox DocuPrint C8	cdj550
Xerox DocuPrint C20	cljet5
Xerox DocuPrint XJ6C	cdj550
Xerox Phaser 3150	pxlmono
Xerox Phaser 8860	cljet5
Xerox WorkCentre 450cp	cdj550
Xerox WorkCentre C2424	cljet5
Xerox WorkCentre PE120	pxlmono

## Brother

Printer	Driver
Brother DCP-1200	ljet2p
Brother DCP-7010	ljet4, ljet4d
Brother DCP-7020	ljet4, ljet4d
Brother DCP-8045D	lj5gray, lj5mono, ljet4, pxlmono
Brother HL-4Ve	laserjet
Brother HL-8	ljetplus
Brother HL-10h	ljet4
Brother HL-10V	ljet3
Brother HL-630	ljet2p
Brother HL-660	ljet4
Brother HL-760	ljet4
Brother HL-960	ljet4
Brother HL-1040	ljet2p
Brother HL-1050	ljet4
Brother HL-1060	ljet4
Brother HL-1070	ljet4
Brother HL-1240	laserjet
Brother HL-1250	lj5gray, lj5mono, ljet4, pxlmono
Brother HL-1260	ljet4
Brother HL-1270N	lj5gray, lj5mono, ljet4, pxlmono
Brother HL-1430	ljet4
Brother HL-1440	ljet4
Brother HL-1450	lj5gray, lj5mono, ljet4, pxlmono
Brother HL-1470N	lj5gray, lj5mono, ljet4, pxlmono
Brother HL-1650	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Brother HL-1660e	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Brother HL-1670N	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Brother HL-1850	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Brother HL-1870N	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Brother HL-2030 (new)	lj5gray, ljet4
Brother HL-2035 (new)	lj5gray, ljet4
Brother HL-2040 (new)	ljet4
Brother HL-2060	lj5gray, lj5mono, ljet4, pxlmono
Brother HL-2140 (new)	lj5gray, ljet4, pxlmono
Brother HL-2170W (new)	ljet4, ljet4d
Brother HL-2400CeN	lj5gray, lj5mono, pxlcolor
Brother HL-2460	lj5gray, lj5mono, pxlmono
Brother HL-2460N	lj5gray, lj5mono, pxlmono
Brother HL-2600CN	lj5gray, lj5mono, pxlcolor
Brother HL-3260N	lj5gray, pxlmono
Brother HL-3400CN	lj5gray, lj5mono, pxlcolor
Brother HL-3450CN	lj5gray, lj5mono, pxlcolor
Brother HL-4000CN	lj5gray, lj5mono, pxlcolor

**Printer list**  
**Brother**

Printer	Driver
Brother HL-4040CN (new)	pxlcolor
Brother HL-4050CDN (new)	pxlcolor
Brother HL-4070CDW (new)	pxlcolor
Brother HL-5030	lj5gray, lj5mono, ljet4, pxlmono
Brother HL-5040	lj5gray, lj5mono, ljet4, pxlmono
Brother HL-5050	lj5gray, lj5mono, ljet4, pxlmono
Brother HL-5070N	lj5gray, lj5mono, ljet4, pxlmono
Brother HL-5140	lj5gray, lj5mono, ljet4, pxlmono
Brother HL-5150D	lj5gray, lj5mono, ljet4, pxlmono
Brother HL-5170DN	lj5gray, lj5mono, ljet4, pxlmono
Brother HL-5250DN	lj5gray, pxlmono
Brother HL-7050	lj5gray, lj5mono, pxlmono
Brother HL-7050N	lj5gray, lj5mono, pxlmono
Brother MC-3000 (new)	epson
Brother MFC-6550MC	ljet2p
Brother MFC-7150C (new)	epsonc
Brother MFC-8300	laserjet, ljet3
Brother MFC-8600	laserjet, ljet3
Brother MFC-9100c (new)	epsonc
Brother MFC-9500	ljet2p
Brother MFC-9600	ljet2p
Brother MFX-9840CDW (new)	pxlcolor

## Canon

Printer	Driver
Canon BJ-5 (new)	bj10e
Canon BJ-10e (new)	bj10e
Canon BJ-20 (new)	bj10e
Canon BJ-30 (new)	bj200
Canon BJ-100 (new)	bj200
Canon BJ-200 (new)	bj200
Canon BJC-210 (new)	bjc600
Canon BJC-240 (new)	bjc600
Canon BJC-250 (new)	bjc600
Canon BJ-330 (new)	bj200
Canon BJC-50 (new)	bjc600
Canon BJC-55 (new)	bjc600
Canon BJC-70 (new)	bj200, bjc600
Canon BJC-80 (new)	bjc600
Canon BJC-85 (new)	bjc600
Canon BJC-210 (new)	bj200
Canon BJC-210SP (new)	bj200
Canon BJC-240 (new)	bj200
Canon BJC-250 (new)	bj200
Canon BJC-600 (new)	bjc600
Canon BJC-610 (new)	bjc600
Canon BJC-620 (new)	bjc600
Canon BJC-800 (new)	bjc800
Canon BJC-1000 (new)	bj200, bjc600
Canon BJC-2000 (new)	bjc600
Canon BJC-2010 (new)	bjc600
Canon BJC-2100 (new)	bjc600
Canon BJC-2110 (new)	bjc600
Canon BJC-3000 (new)	bjc600
Canon BJC-4000 (new)	bjc600
Canon BJC-4100 (new)	bjc600
Canon BJC-4200 (new)	bjc600
Canon BJC-4300 (new)	bjc600
Canon BJC-4310SP (new)	bjc600
Canon BJC-4400 (new)	bjc600
Canon BJC-4550 (new)	bj200, bjc800
Canon BJC-7000 (new)	bjc800
Canon BJC-7100 (new)	bjc800
Canon BJC-4550 (new)	bjc800
Canon GP 335	ljet4
Canon i450 (new)	bjc800
Canon imageRunner 330s	ljet4
Canon imageRunner 1023 (new)	pxlmono

**Printer list**  
Canon

Printer	Driver
Canon imageRunner 1023iF (new)	pxlmono
Canon imageRunner 1023N (new)	pxlmono
Canon imageRunner 2016 (new)	pxlmono
Canon imageRunner 2016i (new)	pxlmono
Canon imageRunner 2018 (new)	pxlmono
Canon imageRunner 2020 (new)	pxlmono
Canon imageRunner 2020i (new)	pxlmono
Canon imageRunner 2022 (new)	pxlmono
Canon imageRunner 2200	ljet4, ljet4d, pxlmono
Canon imageRunner 2230 (new)	pxlmono
Canon imageRunner 2270 (new)	pxlmono
Canon imageRunner 2800	ljet4, ljet4d, pxlmono
Canon imageRunner 3025 (new)	pxlmono
Canon imageRunner 3225 (new)	pxlmono
Canon imageRunner 3300	ljet4, ljet4d, pxlmono
Canon imageRunner 3570	ljet4, ljet4d, pxlmono
Canon imageRunner 4570	ljet4, ljet4d, pxlmono
Canon imageRunner 5000	ljet4, ljet4d, pxlmono
Canon imageRunner 5570	ljet4, ljet4d, pxlmono
Canon imageRunner 6000	ljet4, ljet4d, pxlmono
Canon imageRunner 6570	ljet4, ljet4d, pxlmono
Canon imageRunner 7086 (new)	pxlmono
Canon imageRunner 7095 (new)	pxlmono
Canon imageRunner 7105 (new)	pxlmono
Canon imageRunner 8070 (new)	pxlmono
Canon imageRunner 8500	ljet4, ljet4d
Canon imageRunner C2550 (new)	pxlcolor
Canon imageRunner C3200 (new)	cljet5
Canon imageRunner C5058 (new)	pxlcolor
Canon imageRunner C5068 (new)	pxlcolor
Canon imageRunner C5180 (new)	cljet5
Canon imageRunner C5185 (new)	pxlcolor
Canon imageRunner C5185i (new)	pxlcolor
Canon imageRunner C5870U	ljet4, ljet4d
Canon imageRunner C6870U	ljet4, ljet4d
Canon LBP-4sx	ljet3
Canon LBP-430	ljet3, ljet4
Canon LBP-1000	lj5gray, lj5mono, ljet4, pxlmono
Canon LBP-1260	ljet4
Canon LBP-1760	lj5gray, lj5mono, ljet4, pxlmono
Canon LBP-3360	ljet4, ljet4d, pxlmono
Canon LBP-3460 (new)	pxlmono
Canon LBP-5360 (new)	cljet5
Canon LBP-5960 (new)	cljet5
Canon LBP-5970 (new)	pxlmono

**Printer list**  
Canon

Printer	Driver
Canon LBP-5975 (new)	pxlmono
Canon MultiPASS C2500 (new)	bjc600
Canon MultiPASS C3000 (new)	bjc600
Canon MultiPASS C3500 (new)	bjc600
Canon MultiPASS C5000 (new)	bjc600
Canon MultiPASS C5500 (new)	bjc600
Canon S100 (new)	bjc600
Canon S400 (new)	bjc600

## Epson

Printer	Driver
Epson ActionLaser II	laserjet
Epson ActionLaser 1100	ljet3
Epson AcuLaser C1900	ljet4
Epson AcuLaser C1900 PS3 (new)	ljet4
Epson AcuLaser C2000	lj5gray, lj5mono, ljet4
Epson AcuLaser C2000 PS3 (new)	lj5gray, ljet4
Epson AcuLaser C4000	ljet4
Epson AcuLaser C4000 PS3 (new)	ljet4
Epson AcuLaser C4100	ljet4
Epson AcuLaser C8500	lj5gray, lj5mono, ljet4
Epson AcuLaser C8500PS	lj5gray, lj5mono, ljet4
Epson AcuLaser C8600	lj5gray, lj5mono, ljet4
Epson AcuLaser C8600 PS3 (new)	lj5gray, ljet4, jlet4d
Epson Dot Matrix (new)	epson, epsonc, eps9mid, eps9high
Epson EPL-5200	ljet3
Epson EPL-5200+	ljet3
Epson EPL-5700	ljet4
Epson EPL-5800PS	lj5gray, lj5mono, ljet4, pxlmono
Epson EPL-5900 PS3 (new)	lj5gray, ljet4, pxlmono
Epson EPL-6100	lj5gray, lj5mono, ljet4, pxlmono
Epson EPL-6100 PS3 (new)	lj5gray, ljet4, pxlmono
Epson EPL-7100	laserjet
Epson EPL-N1600	ljet4
Epson EPL-N1600PS	ljet4
Epson EPL-N2050	lj5gray, lj5mono, ljet4, pxlmono
Epson EPL-N2050+	lj5gray, lj5mono, ljet4, pxlmono
Epson EPL-N2050PS	lj5gray, lj5mono, ljet4, pxlmono
Epson EPL-N2050PS+	lj5gray, lj5mono, ljet4, pxlmono
Epson EPL-N2120	lj5gray, lj5mono, ljet4, pxlmono
Epson EPL-N2500	lj5gray, lj5mono, ljet4, pxlmono
Epson EPL-N2500 PS3 (new)	lj5gray, ljet4, pxlmono
Epson EPL-N2750	lj5gray, lj5mono, ljet4, pxlmono
Epson EPL-N2750PS	lj5gray, lj5mono, ljet4, pxlmono
Epson L1000 (new)	epson
Epson LQ24 (new)	epsonc
Epson LQ-500 (new)	epson
Epson LQ-570+ (new)	epson
Epson LQ-2550	epsonc
Epson LX-300+ (USB)	lx-300, epson, ibmpro
Epson LX-1050 (new)	epson

## Fujitsu

Printer	Driver
Fujitsu 1200	epsonc
Fujitsu 2400	epsonc
Fujitsu 3400	epsonc
Fujitsu PrintPartner 10V	laserjet, ljet4
Fujitsu PrintPartner 16DV	ljet4
Fujitsu PrintPartner 20W	lj5gray, lj5mono, ljet4, pxlmono
Fujitsu PrintPartner 8000	ljet3

## Generic

Printer	Driver
Generic ESC/P Dot Matrix Printer (new)	epson, epsonc, eps9mid, eps9high
Generic IMB-Compatible Dot Matrix Printer (new)	ibmpro
Generic PCL 4 Printer	laserjet, ljet2p, ljetplus
Generic PCL 5 Printer	ljet3, ljet3d
Generic PCL 5c Printer	ljet4, ljet4d, cljet5
Generic PCL 5e Printer	ljet4, ljet4d
Generic PCL 6/PCL XL Printer	lj5gray, lj5mono, ljet4, ljet4d, pxlcolor, pxlmono

# HP

Printer	Driver
HP - 1005	ljet4
HP 2550	ljet4
HP Business Inkjet 3000 (new)	pxlcolor
HP Color LaserJet 5 (new)	cljet5
HP Color LaserJet 5M (new)	cljet5
HP Color LaserJet 4550	lj5gray, lj5mono, ljet4
HP Color LaserJet 4600	lj5gray, lj5mono, ljet4
HP Color LaserJet 5000	ljet4, cljet5
HP Color LaserJet 5500 (new)	pxlcolor
HP DeskJet	deskjet
HP DeskJet 200	deskjet
HP Deskjet 310 (new)	cdj550
HP DeskJet 320	djet500, cdj550
HP DeskJet 340C	deskjet
HP Deskjet f340 (new)	deskjet
HP DeskJet 500	djet500
HP DeskJet 510	djet500
HP DeskJet 520	djet500
HP DeskJet 540	djet500
HP DeskJet 540C	djet500, cdj550
HP DeskJet 550C	djet500, cdj550
HP Deskjet 560C (new)	cdj550
HP DeskJet 600	djet500, cdj550
HP Deskjet 610CL (new)	cdj550
HP Deskjet 612C (new)	cdj550
HP Deskjet 648C (new)	cdj550
HP Deskjet 692C (new)	cdj550
HP Deskjet 694C (new)	cdj550
HP Deskjet 812C (new)	cdj550
HP Deskjet 840C (new)	cdj550
HP Deskjet 842C (new)	cdj550
HP Deskjet 930C (new)	cdj550
HP Deskjet 932C (new)	cdj550
HP Deskjet 950C (new)	cdj550
HP Deskjet 970C (new)	cdj550
HP Deskjet 990C (new)	cdj550
HP Deskjet 1200C (new)	pjl300
HP Deskjet 2566	djet500
HP Deskjet D2466	djet500
HP Deskjet D2563	djet500
HP DeskJet Plus	deskjet
HP Deskjet Series	ljet4
HP Laser M1120 MFP	ljet4

Printer	Driver
HP LaserJet	laserjet
HP LaserJet Plus	ljetplus
HP LaserJet 2	laserjet
HP LaserJet 2D	ljet2p
HP LaserJet 2P	ljet2p
HP LaserJet 2P Plus	ljet2p
HP LaserJet 3	ljet3
HP LaserJet 3D	ljet3d
HP LaserJet 3P w/PS	ljet3, ljet3d
HP LaserJet 3P w/ PCL5	ljet4
HP LaserJet 4	ljet4
HP LaserJet 4 Plus	ljet4
HP LaserJet 4L	ljet4
HP LaserJet 4M	ljet4, ljet4d
HP LaserJet 4ML	ljet4
HP LaserJet 4P	ljet4
HP LaserJet 4Si	ljet4, ljet4d
HP LaserJet 4V	ljet4
HP LaserJet 5	lj5gray, lj5mono, ljet4, pxlmono
HP LaserJet 5L	ljet4
HP LaserJet 5P	ljet4
HP LaserJet 5Si	lj5gray, lj5mono, ljet4
HP LaserJet 6	lj5gray, lj5mono, ljet4
HP LaserJet 6L	ljet4
HP LaserJet 6P	ljet4
HP LaserJet 1015	ljet4
HP LaserJet 1100	ljet4
HP LaserJet 1100A	ljet4
HP LaserJet 1150	ljet4
HP LaserJet 1200	lj5gray, lj5mono, ljet4, pxlmono
HP LaserJet 1220	lj5gray, lj5mono, ljet4, pxlmono
HP LaserJet 1300	lj5gray, lj5mono, ljet4, pxlmono
HP LaserJet 1320 (new)	pxlmono
HP LaserJet 2100	lj5gray, lj5mono, ljet4, pxlmono
HP LaserJet 2100M	lj5gray, lj5mono, ljet4, pxlmono
HP LaserJet 2200	lj5gray, lj5mono, ljet4, pxlmono
HP LaserJet 2300	lj5gray, lj5mono, ljet4, pxlmono
HP LaserJet 2500 (new)	pxlcolor
HP LaserJet 3200	lj5gray, lj5mono, ljet4, pxlmono
HP LaserJet 3200m	lj5gray, lj5mono, ljet4, pxlmono
HP LaserJet 3200se	lj5gray, lj5mono, ljet4, pxlmono
HP LaserJet 3300 MFP	lj5gray, lj5mono, pxlmono
HP LaserJet 3310 MFP	lj5gray, lj5mono, pxlmono
HP LaserJet 3320 MFP	lj5gray, lj5mono, pxlmono
HP LaserJet 3320N MFP	lj5gray, lj5mono, pxlmono

**Printer list**  
HP

Printer	Driver
HP LaserJet 3330 MFP	lj5gray, lj5mono, pxlmono
HP LaserJet 3390 (new)	lj5gray, ljet4, ljet4d, pxlmono
HP LaserJet 3392 (new)	lj5gray, ljet4, ljet4d, pxlmono
HP LaserJet 4000	lj5gray, lj5mono, pxlmono
HP LaserJet 4050	lj5gray, lj5mono, ljet4, pxlmono
HP LaserJet 4100	lj5gray, lj5mono, ljet4, pxlmono
HP LaserJet 4200 (new)	pxlmono
HP LaserJet 4300 (new)	pxlmono
HP LaserJet 4350	lj5gray, pxlmono
HP LaserJet 4500 (new)	pxlcolor
HP LaserJet 4550 (new)	pxlcolor
HP LaserJet 4600 (new)	pxlcolor
HP LaserJet 5000 (new)	pxlcolor
HP LaserJet 5100 (new)	pxlmono
HP LaserJet 5500 (new)	pxlcolor
HP LaserJet 8150	lj5gray, lj5mono, ljet4, pxlmono
HP LaserJet 9000	lj5gray, lj5mono, ljet4, pxlmono
HP LaserJet M1522nf MFP (new)	pxlmono
HP LaserJet M1522n MFP(new)	pxlmono
HP new models	ljet4
HP OfficeJet	djet500
HP OfficeJet R45	deskjet
HP OfficeJet 300	djet500
HP OfficeJet 330	djet500
HP OfficeJet 350	djet500
HP OfficeJet 500 (new)	cdj550
HP OfficeJet 600 (new)	cdj550
HP OfficeJet G55 (new)	cdj550
HP OfficeJet G85 (new)	cdj550
HP OfficeJet G95 (new)	cdj550
HP OfficeJet LX	djet500
HP OfficeJet Pro 1150C (new)	cdj550
HP OfficeJet Pro 1170C (new)	cdj550
HP OfficeJet R45 (new)	cdj550
HP OfficeJet T45 (new)	cdj550
HP OfficeJet T65 (new)	cdj550
HP - P2035	ljet4
HP PaintJet (new)	pj
HP PaintJet XL (new)	pjxl
HP PaintJet XL300 (new)	pjxl300
HP PhotoSmart P100 (new)	pjxl300
HP PhotoSmart P1000 (new)	cdj550
HP PhotoSmart P1100 (new)	cdj550

## IBM

Printer	Driver
IBM Execjet 4072 (new)	bj200
IBM Infoprint 12	lj5gray, lj5mono, ljet4, pxlmono
IBM Page Printer 3112	ljet4
IBM ProPrinterII (new)	ibmpro
IBM 4019	laserjet
IBM 4029 030 LaserPrinter 10	laserjet, ljet3
IBM 4312 (new)	laserjet

## Kyocera

Printer	Driver
Kyocera Ci-1100 (new)	cljet5
Kyocera CS-250ci KPD L (new)	cljet5
Kyocera CS-300ci KPD L (new)	cljet5
Kyocera CS-400ci KPD L (new)	cljet5
Kyocera CS-500ci KPD L (new)	cljet5
Kyocera CS-1815	ljet4, ljet4d
Kyocera CS-C2525E KPD L (new)	cljet5
Kyocera CS-C3225E KPD L (new)	cljet5
Kyocera CS-C3232E KPD L (new)	cljet5
Kyocera CS-C4035E KPD L (new)	cljet5
Kyocera F-800	laserjet, ljetplus
Kyocera F-800T	laserjet, ljetplus
Kyocera F-820	laserjet, ljetplus
Kyocera F-1000	laserjet, ljetplus
Kyocera F-1010	laserjet, ljetplus
Kyocera F-1200	laserjet, ljetplus
Kyocera F-1200S	laserjet, ljetplus
Kyocera F-1800	laserjet, ljetplus
Kyocera F-2000	laserjet, ljetplus
Kyocera F-2010	laserjet, ljetplus
Kyocera F-2200	laserjet, ljetplus
Kyocera F-2200S	laserjet, ljetplus
Kyocera F-3000	laserjet, ljetplus
Kyocera F-3010	laserjet, ljetplus
Kyocera F-3300	laserjet, ljetplus
Kyocera F-5000	laserjet, ljetplus
Kyocera FS-400	laserjet, ljet3
Kyocera FS-600	ljet4, ljet4d
Kyocera FS-600 (KPD L-2)	ljet4, ljet4d
Kyocera FS-680	ljet4, ljet4d
Kyocera FS-800	ljet4, ljet4d
Kyocera FS-850	laserjet, ljet3
Kyocera FS-920	ljet4, ljet4d
Kyocera FS-1000	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-1000+	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-1010	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-1018MFP	ljet4, ljet4d
Kyocera FS-1020D	ljet4, ljet4d
Kyocera FS-1030D	ljet4, ljet4d
Kyocera FS-1050	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-1118 MFP	ljet4, ljet4d
Kyocera FS-1200	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-1500	laserjet, ljet3, ljet3d

**Printer list**  
Kyocera

Printer	Driver
Kyocera FS-1550	laserjet, ljet3, ljet3d
Kyocera FS-1550+	laserjet, ljet3, ljet3d
Kyocera FS-1600	ljet4, ljet4d
Kyocera FS-1600+	ljet4, ljet4d
Kyocera FS-1700	ljet4, ljet4d
Kyocera FS-1700+	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-1714M	ljet4, ljet4d
Kyocera FS-1750	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-1800	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-1800+	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-1900	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-1920	ljet4, ljet4d
Kyocera FS-2000D	ljet4, ljet4d
Kyocera FS-3400	laserjet, ljet3, ljet3d
Kyocera FS-3400+	laserjet, ljet3, ljet3d
Kyocera FS-3500	laserjet, ljet3, ljet3d
Kyocera FS-3600	ljet4, ljet4d
Kyocera FS-3600+	ljet4, ljet4d
Kyocera FS-3700	ljet4, ljet4d
Kyocera FS-3700+	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-3718M	ljet4, ljet4d
Kyocera FS-3750	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-3800	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-3820N	ljet4, ljet4d
Kyocera FS-3830N	ljet4, ljet4d
Kyocera FS-3900DN	ljet4, ljet4d
Kyocera FS-4000DN	ljet4, ljet4d
Kyocera FS-5500	laserjet, ljet3
Kyocera FS-5800C	ljet4, ljet4d, cljet5
Kyocera FS-5900C	ljet4, ljet4d, cljet5
Kyocera FS-6020	ljet4, ljet4d
Kyocera FS-6026	ljet4, ljet4d
Kyocera FS-6300	ljet4, ljet4d
Kyocera FS-6500	ljet4, ljet4d
Kyocera FS-6500+	ljet4, ljet4d
Kyocera FS-6700	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-6750	ljet4, ljet4d
Kyocera FS-6900	ljet4, ljet4d
Kyocera FS-6950DN	ljet4, ljet4d
Kyocera FS-7000	ljet4, ljet4d
Kyocera FS-7000+	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-7028M	ljet4, ljet4d
Kyocera FS-8000C	ljet4, ljet4d, cljet5
Kyocera FS-9000	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-9100DN	lj5gray, lj5mono, ljet4, ljet4d, pxlmono

Printer	Driver
Kyocera FS-9130DN	ljet4, ljet4d
Kyocera FS-9500DN	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera FS-9530DN	ljet4, ljet4d
Kyocera FS-C5015N (new)	cljet5
Kyocera FS-C5016N (new)	cljet5
Kyocera FS-C5020N (new)	cljet5
Kyocera FS-C5025N (new)	cljet5
Kyocera FS-C5030N (new)	cljet5
Kyocera FS-C5100DN KPD (new)	cljet5
Kyocera FS-C5200DN KPD (new)	cljet5
Kyocera FS-C5300DN KPD (new)	cljet5
Kyocera FS-C5350DN KPD (new)	cljet5
Kyocera FS-C5400DN KPD (new)	cljet5
Kyocera FS-C8008N (new)	cljet5
Kyocera FS-C8100DN (new)	cljet5
Kyocera FS-C8100DN+KPD (new)	cljet5
Kyocera KM-1510	ljet4, ljet4d
Kyocera KM-1530	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera KM-1810	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera KM-1815	ljet4, ljet4d
Kyocera KM-1820	ljet4, ljet4d
Kyocera KM-2030	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera KM-2530	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera KM-3050	ljet4, ljet4d
Kyocera KM-3530	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera KM-4050	ljet4, ljet4d
Kyocera KM-4230	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera KM-4230/5230	ljet4, ljet4d
Kyocera KM-4530	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera KM-5050	ljet4, ljet4d
Kyocera KM-5230	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera KM-5530	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera KM-6030	ljet4, ljet4d
Kyocera KM-6230	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Kyocera KM-8030	ljet4, ljet4d
Kyocera KM-C830 (new)	cljet5
Kyocera KM-C850 (new)	cljet5
Kyocera KM-C2520 (new)	cljet5
Kyocera KM-C2525E KPD (new)	cljet5
Kyocera KM-C2630 (new)	cljet5
Kyocera KM-C3225 (new)	cljet5
Kyocera KM-C3225E KPD (new)	cljet5
Kyocera KM-C3232 (new)	cljet5
Kyocera KM-C3232E KPD (new)	cljet5
Kyocera KM-C4035E KPD (new)	cljet5

**Printer list**  
Kyocera

Printer	Driver
Kyocera P-2000	laserjet, ljetplus
Kyocera P-2002	laserjet, ljetplus
Kyocera TASKalfa 250ci KPD (new)	cljet5
Kyocera TASKalfa 300ci KPD (new)	cljet5
Kyocera TASKalfa 400ci KPD (new)	cljet5
Kyocera TASKalfa 500ci KPD (new)	cljet5

## Lexmark

Printer	Driver
Lexmark 3000 (new)	cdj550
Lexmark C720n (new)	pxlcolor
Lexmark E120	ljet4, ljet4d
Lexmark E120n (new)	ljet4, ljet4d, pxlmono
Lexmark E230 (new)	pxlmono
Lexmark E238 (new)	pxlmono
Lexmark E320	ljet4, ljet4d
Lexmark E350d	lj5gray, pxlmono
Lexmark E352dn	lj5gray, pxlmono
Lexmark E450dn (new)	lj5gray, pxlmono
Lexmark Optra E	ljet4
Lexmark Optra E+	ljet4
Lexmark Optra E321	ljet4
Lexmark Valuewriter 300	ljet2p
Lexmark Winwriter 400	ljetplus
Lexmark X342n (new)	pxlmono
Lexmark X642e (new)	pxlmono

## Minolta

Printer	Driver
Minolta PagePro 6	ljet4
Minolta PagePro 6e	ljet4
Minolta PagePro 6ex	ljet4
Minolta PagePro 8	ljet4
Minolta PagePro 8L	ljet2p
Minolta PagePro 1100	ljet4, pxlmono

## NEC

Printer	Driver
NEC PinWriter P6 (new)	epsonc
NEC PinWriter P6 plus (new)	epsonc
NEC SuperScript 660i	ljet4
NEC SuperScript 860	ljet2p
NEC SuperScript 870	ljet2p
NEC SuperScript 1260	ljet2p
NEC SuperScript 1400	ljet4
NEC SuperScript 1800	ljet4

## Oce

Printer	Driver
Oce 8445	ljet3, ljet3d
Oce 8465	ljet3, ljet3d
Oce fx2080 (new)	pxlmono
Oce sx 1480 (new)	pxlmono
Oce VarioPrint 2050	ljet3, ljet3d
Oce VarioPrint 2060	ljet3, ljet3d
Oce VarioPrint 2070	ljet3, ljet3d
Oce VarioPrint 2090	ljet3, ljet3d
Oce VarioPrint2105	ljet3, ljet3d

## Okidata

Printer	Driver
OKI B410 (new)	ljet4, ljet4d, pxlmono
OKI B2400 (new)	pxlmono
OKI B2400n (new)	pxlmono
OKI B4200 (new)	pxlmono
OKI B4250 (new)	pxlmono
OKI B4300	lj5gray, pxlmono
OKI B4350 (new)	ljet4, pxlmono
OKI B4400 (new)	pxlmono
OKI B4400n (new)	pxlmono
OKI B4600 (new)	pxlmono
OKI B4600n (new)	pxlmono
OKI B4600n PS (new)	pxlmono
OKI B4600 PS (new)	pxlmono
OKI B6300 (new)	pxlmono
OKI C5300 (new)	cljet5
OKI C5400 (new)	cljet5
OKI C6000 (new)	pxlcolor
OKI C6100 (new)	pxlcolor
OKI C7300 (new)	cljet5
OKI C7350 (new)	pxlcolor
OKI C7550 (new)	pxlcolor
OKI C8800 (new)	pxlcolor
OKI C9200 (new)	cljet5
OKI C9400 (new)	cljet5
OKI C9650 (new)	pxlcolor
OKI C9800 (new)	pxlcolor
OKI Microline 192+ (new)	ibmpro
OKI Microline 320 turbo	okiibm
OKI Microline IBM compatible 9 pin (new)	okiibm
OKI ML 320 (new)	okiibm
OKI ML 321 (new)	okiibm
OKI ML 380 (new)	epson
OKI Okijet 2500 (new)	cdj550
OKI Okipage 6e	ljet4
OKI Okipage 6ex	ljet4
OKI Okipage 8p	ljet4
OKI Okipage 10e	ljet4
OKI Okipage 10ex	lj5gray, lj5mono, ljet4, pxlmono
OKI Okipage 14ex	lj5gray, lj5mono, ljet4, pxlmono
OKI OL400	laserjet
OKI OL400e	ljet2p
OKI OL400ex	ljet2p
OKI OL410e	ljet4

**Printer list**  
Okidata

Printer	Driver
OKI OL600e	ljet2p
OKI OL610e/S	laserjet
OKI OL800	laserjet
OKI OL810ex	ljet4
OKI Super 6e	ljet4

## Panasonic

Printer	Driver
Panasonic DP-180 (new)	pxlcolor
Panasonic DP-2330 (new)	pxlmono
Panasonic DP-4510 (new)	pxlmono
Panasonic DP-4520 (new)	pxlmono
Panasonic DP-8016P (new)	pxlmono
Panasonic DP-8020E (new)	pxlmono
Panasonic DP-8045 (new)	pxlmono
Panasonic DP-8060 (new)	pxlmono
Panasonic DP-C213 (new)	pxlmono
Panasonic DP-C262 (new)	pxlmono
Panasonic DP-C264 (new)	pxlmono
Panasonic DP-C265 (new)	pxlcolor
Panasonic DP-C305 (new)	pxlcolor
Panasonic DP-C322 (new)	pxlcolor
Panasonic DP-C354 (new)	pxlmono
Panasonic DP-C405 (new)	pxlcolor
Panasonic KX-MC6020 (new)	pxlcolor
Panasonic KX-P1123 (new)	epson
Panasonic KX-P1124 (new)	epson
Panasonic KX-P1150 (new)	eps9high
Panasonic KX-P1624 (new)	epson
Panasonic KX-P2023 (new)	epson
Panasonic KX-P2123 (new)	epson
Panasonic KX-P2135 (new)	epsonc
Panasonic KX-P2150 (new)	ibmpro
Panasonic KX-P4410	ljet2p
Panasonic KX-P4450	ljet3
Panasonic KX-P6150	laserjet
Panasonic KX-P6500	ljet2p

## Samsung

Printer	Driver
Samsung CLP-325w (new)	pxlcolor
Samsung CLP-350 (new)	pxlcolor
Samsung CLP-550 (new)	pxlcolor
Samsung CLP-550N (new)	pxlcolor
Samsung CLP-650 (new)	pxlcolor
Samsung CLP-650N (new)	pxlcolor
Samsung CLP-660N (new)	pxlcolor
Samsung CLP-660ND (new)	pxlcolor
Samsung ML-85	ljet4
Samsung ML-1250	lj5gray, lj5mono, ljet4, pxlmono
Samsung ML-1450	ljet4, pxlmono
Samsung ML-1450PS	ljet4, pxlmono
Samsung ML-1650	lj5gray, lj5mono, ljet4, pxlmono
Samsung ML-1651N	lj5gray, lj5mono, ljet4, pxlmono
Samsung ML-1750	ljet4, pxlmono
Samsung ML-2150	ljet4, ljet4d, pxlmono
Samsung ML-2150PS	ljet4, ljet4d, pxlmono
Samsung ML-2151N	ljet4, ljet4d, pxlmono
Samsung ML-2151NPS	ljet4, ljet4d, pxlmono
Samsung ML-2152W	ljet4, ljet4d, pxlmono
Samsung ML-2152WPS	ljet4, ljet4d, pxlmono
Samsung ML-2250	ljet4, ljet4d, pxlmono
Samsung ML-2550	ljet4, ljet4d, pxlmono
Samsung ML-2551N	ljet4, ljet4d, pxlmono
Samsung ML-2552W	ljet4, ljet4d, pxlmono
Samsung ML-2570 (new)	pxlmono
Samsung ML-2571N (new)	pxlmono
Samsung ML-2850D (new)	pxlmono
Samsung ML-2851ND (new)	pxlmono
Samsung ML-3051N (new)	pxlmono
Samsung ML-3051ND (new)	pxlmono
Samsung ML-3470N (new)	pxlmono
Samsung ML-3471ND (new)	pxlmono
Samsung ML-4050N (new)	pxlmono
Samsung ML-4550 (new)	pxlmono
Samsung ML-4551N (new)	pxlmono
Samsung ML-4551ND (new)	pxlmono
Samsung ML-4600	lj5gray, lj5mono, ljet4, pxlmono
Samsung ML-5000a	ljet4
Samsung ML-6000	ljet4
Samsung ML-6060 (new)	pxlmono
Samsung ML-6060N (new)	pxlmono
Samsung ML-6060S (new)	pxlmono

**Printer list**  
Samsung

Printer	Driver
Samsung ML-6100	ljet4
Samsung ML-7000	lj5gray, lj5mono, ljet4, pxlmono
Samsung ML-7000N	lj5gray, lj5mono, ljet4, pxlmono
Samsung ML-7000P	lj5gray, lj5mono, ljet4, pxlmono
Samsung ML-7050	lj5gray, lj5mono, ljet4, pxlmono
Samsung ML-7300	ljet4, ljet4d, pxlmono
Samsung ML-7300N	ljet4, ljet4d, pxlmono
Samsung QL-5100A	ljet4
Samsung QL-6050	ljet4

## Sharp

Printer	Driver
Sharp AR-155FG (new)	pxlmono
Sharp AR-160M (new)	pxlmono
Sharp AR-161	lj5gray, lj5mono, ljet4, pxlmono
Sharp AR-163FG (new)	pxlmono
Sharp AR-163G (new)	pxlmono
Sharp AR-168D (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-168S (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-200M (new)	pxlmono
Sharp AR-205FG (new)	pxlmono
Sharp AR-205G (new)	pxlmono
Sharp AR-266FP (new)	pxlmono
Sharp AR-266S (new)	pxlmono
Sharp AR-311FP (new)	pxlmono
Sharp AR-311N (new)	pxlmono
Sharp AR-311S (new)	pxlmono
Sharp AR-351FP (new)	pxlmono
Sharp AR-351N (new)	pxlmono
Sharp AR-351S (new)	pxlmono
Sharp AR-451FP (new)	pxlmono
Sharp AR-451N (new)	pxlmono
Sharp AR-451S (new)	pxlmono
Sharp AR-555M (new)	pxlmono
Sharp AR-555S (new)	pxlmono
Sharp AR-625M (new)	pxlmono
Sharp AR-625S (new)	pxlmono
Sharp AR-705M (new)	pxlmono
Sharp AR-705S (new)	pxlmono
Sharp AR-5220 (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-BC260 (new)	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-BC320 (new)	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C170FP (new)	pxlcolor
Sharp AR-C170M (new)	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C172FP (new)	pxlcolor
Sharp AR-C172M (new)	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C260 (new)	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C260F (new)	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C260FP (new)	pxlcolor
Sharp AR-C260M (new)	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C260P (new)	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C260S (new)	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C261F (new)	ljet4, ljet4d, cljet5, pxlcolor
Sharp AR-C261FP (new)	pxlcolor
Sharp AR-C261M (new)	ljet4, ljet4d, cljet5, pxlcolor

**Printer list**  
Sharp

Printer	Driver
Sharp AR-C261S (new)	ljet4, ljet4d, cljet5, pxcolor
Sharp AR-C262FP (new)	pxlcolor
Sharp AR-C262M (new)	ljet4, ljet4d, cljet5, pxcolor
Sharp AR-C262S (new)	ljet4, ljet4d, cljet5, pxcolor
Sharp AR-M150 (new)	pxlmono
Sharp AR-M155 (new)	pxlmono
Sharp AR-M160 (new)	pxlmono
Sharp AR-M161 (new)	pxlmono
Sharp AR-M162 (new)	pxlmono
Sharp AR-M165 (new)	pxlmono
Sharp AR-M205 (new)	pxlmono
Sharp AR-M206 (new)	pxlmono
Sharp AR-M207 (new)	pxlmono
Sharp AR-M236 (new)	pxlmono
Sharp AR-M237 (new)	pxlmono
Sharp AR-M276 (new)	pxlmono
Sharp AR-M277 (new)	pxlmono
Sharp AR-M351N (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M351U (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M355N (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M355U (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M451N (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M451U (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M455N (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M455U (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M550N (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M550U (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M620N (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M620U (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M700N (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp AR-M700U (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp ARP350 (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp ARP450 (new)	lj5gray, ljet4, ljet4d, pxlmono
Sharp MX-2300FG (new)	pxlcolor
Sharp MX-2300G (new)	pxlcolor
Sharp MX-2300N (new)	pxlcolor
Sharp MX-2700FG (new)	pxlcolor
Sharp MX-2700G (new)	pxlcolor
Sharp MX-2700N (new)	pxlcolor
Sharp MX-3500FN (new)	pxlcolor
Sharp MX-3500N (new)	pxlcolor
Sharp MX-3501FN (new)	pxlcolor
Sharp MX-3501N (new)	pxlcolor
Sharp MX-4500FN (new)	pxlcolor
Sharp MX-4500N (new)	pxlcolor

**Printer list**  
Sharp

Printer	Driver
Sharp MX-4501FN (new)	pxlcolor
Sharp MX-4501N (new)	pxlcolor

## Star

Printer	Driver
Star LC24-200 (new)	epson, epsonc
Star LC 90 (new)	epson
Star LS-04	ljet2p
Star LaserPrinter 8	laserjet
Star NL-10 (new)	epson
Star NX-1001 (new)	eps9mid, eps9high

## Toshiba

Printer	Driver
Toshiba e-Studio 162 (new)	pxlmono
Toshiba e-Studio 162d (new)	pxlmono
Toshiba e-Studio 167 (new)	pxlmono
Toshiba e-Studio 203l (new)	pxlmono
Toshiba e-Studio 220cp (new)	pxlcolor
Toshiba e-Studio 233 (new)	pxlmono
Toshiba e-Studio 237 (new)	pxlmono
Toshiba e-Studio 283 (new)	pxlmono
Toshiba e-Studio 353 (new)	pxlmono
Toshiba e-Studio 450s (new)	pxlmono
Toshiba e-Studio 453 (new)	pxlmono
Toshiba e-Studio 500p (new)	pxlmono
Toshiba e-Studio 500s (new)	pxlmono
Toshiba e-Studio 523 (new)	pxlmono
Toshiba e-Studio 523t (new)	pxlmono
Toshiba e-Studio 603 (new)	pxlmono
Toshiba e-Studio 603t (new)	pxlmono
Toshiba e-Studio 723 (new)	pxlmono
Toshiba e-Studio 723t (new)	pxlmono
Toshiba e-Studio 853 (new)	pxlmono
Toshiba e-Studio 901 (new)	pxlmono
Toshiba e-Studio 1101 (new)	pxlmono
Toshiba e-Studio 1351 (new)	pxlmono
Toshiba e-Studio 2830c (new)	pxlcolor
Toshiba e-Studio 3530c (new)	pxlcolor
Toshiba e-Studio 4500c (new)	pxlcolor
Toshiba e-Studio 4520c (new)	pxlcolor
Toshiba e-Studio 5500c (new)	pxlcolor

## Xerox

Printer	Driver
Xerox Able 1406	ljet4
Xerox Document Centre 400	ljet4
Xerox Document Homecentre (new)	cdj550
Xerox DocuPrint 4508	ljet4
Xerox DocuPrint C6 (new)	cdj550
Xerox DocuPrint C8 (new)	cdj550
Xerox DocuPrint C20 (new)	cljet5
Xerox DocuPrint N4512	lj5gray, lj5mono, ljet4
Xerox DocuPrint N4512PS	lj5gray, lj5mono, ljet4
Xerox DocuPrint P8e	ljet4
Xerox DocuPrint P12	ljet2p
Xerox DocuPrint P1202	lj5gray, lj5mono, ljet4
Xerox DocuPrint XJ6C (new)	cdj550
Xerox Phaser 2135	ljet4, ljet4d
Xerox Phaser 3150 (new)	pxlmono
Xerox Phaser 3500	lj5gray, pxlmono
Xerox Phaser 4400B	ljet4, ljet4d
Xerox Phaser 4400DT	ljet4, ljet4d
Xerox Phaser 4400DX	ljet4, ljet4d
Xerox Phaser 4400N	ljet4, ljet4d
Xerox Phaser 4500B	ljet4, ljet4d
Xerox Phaser 4500DT	ljet4, ljet4d
Xerox Phaser 4500DX	ljet4, ljet4d
Xerox Phaser 4500N	ljet4, ljet4d
Xerox Phaser 4510B	ljet4, ljet4d
Xerox Phaser 4510DT	ljet4, ljet4d
Xerox Phaser 4510DX	ljet4, ljet4d
Xerox Phaser 4510N	ljet4, ljet4d
Xerox Phaser 5500B	ljet4, ljet4d
Xerox Phaser 5500DN	ljet4, ljet4d
Xerox Phaser 5500DT	ljet4, ljet4d
Xerox Phaser 5500DX	ljet4, ljet4d
Xerox Phaser 5500N	ljet4, ljet4d
Xerox Phaser 6130N	ljet4, ljet4d
Xerox Phaser 6180DN	ljet4, ljet4d
Xerox Phaser 6180MFP-D	ljet4, ljet4d
Xerox Phaser 6200B	ljet4, ljet4d
Xerox Phaser 6200DP	ljet4, ljet4d
Xerox Phaser 6200DX	ljet4, ljet4d
Xerox Phaser 6200N	ljet4, ljet4d
Xerox Phaser 6250B	ljet4, ljet4d
Xerox Phaser 6250DP	ljet4, ljet4d
Xerox Phaser 6250DT	ljet4, ljet4d

Printer	Driver
Xerox Phaser 6250DX	ljet4, ljet4d
Xerox Phaser 6250N	ljet4, ljet4d
Xerox Phaser 6300DN	ljet4, ljet4d
Xerox Phaser 6300N	ljet4, ljet4d
Xerox Phaser 6350DP	ljet4, ljet4d
Xerox Phaser 6350DT	ljet4, ljet4d
Xerox Phaser 6350DX	ljet4, ljet4d
Xerox Phaser 6360DN	ljet4, ljet4d
Xerox Phaser 6360DX	ljet4, ljet4d
Xerox Phaser 7300B	ljet4, ljet4d
Xerox Phaser 7300DN	ljet4, ljet4d
Xerox Phaser 7300DT	ljet4, ljet4d
Xerox Phaser 7300DX	ljet4, ljet4d
Xerox Phaser 7300N	ljet4, ljet4d
Xerox Phaser 7400DN	ljet4, ljet4d
Xerox Phaser 7400DT	ljet4, ljet4d
Xerox Phaser 7400DX	ljet4, ljet4d
Xerox Phaser 7400DXF	ljet4, ljet4d
Xerox Phaser 7400N	ljet4, ljet4d
Xerox Phaser 7700DN	ljet4, ljet4d
Xerox Phaser 7700DX	ljet4, ljet4d
Xerox Phaser 7700GX	ljet4, ljet4d
Xerox Phaser 7750B	ljet4, ljet4d
Xerox Phaser 7750DN	ljet4, ljet4d
Xerox Phaser 7750DXF	ljet4, ljet4d
Xerox Phaser 7750GX	ljet4, ljet4d
Xerox Phaser 7760DN	ljet4, ljet4d
Xerox Phaser 7760DX	ljet4, ljet4d
Xerox Phaser 7760GX	ljet4, ljet4d
Xerox Phaser 8400B	ljet4, ljet4d
Xerox Phaser 8400BD	ljet4, ljet4d
Xerox Phaser 8400DP	ljet4, ljet4d
Xerox Phaser 8400DX	ljet4, ljet4d
Xerox Phaser 8400N	ljet4, ljet4d
Xerox Phaser 8500DN	ljet4, ljet4d
Xerox Phaser 8500N	ljet4, ljet4d
Xerox Phaser 8550DP	ljet4, ljet4d
Xerox Phaser 8550DT	ljet4, ljet4d
Xerox Phaser 8550DX	ljet4, ljet4d
Xerox Phaser 8560DN	ljet4, ljet4d
Xerox Phaser 8860 (new)	cljet5
Xerox WorkCentre 450cp (new)	cdj550
Xerox WorkCentre 7345	ljet4, ljet4d, pxlcolor
Xerox WorkCentre C2424 (new)	cljet5
Xerox WorkCentre PE120 (new)	pxlmono

## Others

Printer	Driver
Anitech M24 (new)	epson
Apollo P-2200 (new)	cdj550
Apollo P-2250 (new)	cdj550
Apple LaserWriter Select 360	ljet2p
Citizen ProJet II	laserjet
DEC 1800	ljet3
DEC DECWriter 110i	djet500
DEC DECwriter 500i	djet500
DEC LN17	ljet4
Dell A920 (new)	deskjet
Dell 3100cn (new)	pxlcolor
Deskjet 3050 all-in-one J610a series (new)	deskjet
Infotec 4651 (new)	pxlmono
Infotec 4651 MF	lj5gray, lj5mono, ljet4
Olivetti JP350S	laserjet
Olivetti PG 306	laserjet
Olivetti JP450	djet500
PCPI 1030	laserjet
QMS LPK-100 (new)	eps9mid
Raven LP-410	ljet2p
Ricoh Aficio 220	lj5gray, lj5mono, ljet4, ljet4d, pxlmono
Ricoh Aficio 401	ljet4
Ricoh Aficio 700	ljet4
Seiko SpeedJET 200	laserjet
Tally MT908	ljet3
Tektronix Phaser 750DP	ljet4, ljet4d
Tektronix Phaser 750DX	ljet4, ljet4d
Tektronix Phaser 750N	ljet4, ljet4d
Tektronix Phaser 750P	ljet4, ljet4d