



HORIBA ABX SAS
 Parc Euromédecine - Rue du Caducée
 B.P. 7290
 34184 Montpellier Cedex 4
 FRANCE
 www.horiba.com



LOT L: 17788

REF: A11A01928

04710819C



2016/01/31

(FR) Constituant (EN) Constituent (DE) Bestandteil (ES) Componente (IT) Componente (PT) Componente (DA) Bestanddel (SV) Konstituent (EL) συστατικό (PL) Składnik (NO) Bestanddel	ABX Pentra 400	Pentra C400	Pentra C200	Unité Unit Einheit Unidad Unità Unidade Enhed Enhet Μονάδα Jednostka Enhet	Valeur titrée Assigned value Deklariertes Wert Valor asignado Valori assegnati Valor determinado Tildelt værdi Tildelt värde Αναγραφόμενη τιμή Wartość przypisana Tildelt verd	Intervalle de confiance Confidence range Vertrauensintervalle Intervalos de confianza Livelli di fiducia Intervalo de confiança Konfidensområde Konfidensintervall Εύρος εμπιστοσύνης Zakres ufności Konfidensspekter	
C				U	\oplus	→ ±3 SDI ←	1 SD
IgA ABX Pentra IgA CP (Ref. A11A01923)	•			g/L mg/dL	1.37 137	1.04 – 1.70 104 - 170	0.110 11.0
		•		g/L mg/dL	1.37 137	1.04 – 1.70 104 - 170	0.110 11.0
			•	g/L mg/dL	1.37 137	1.04 – 1.70 104 - 170	0.110 11.0
				g/L mg/dL	1.37 137	1.04 – 1.70 104 - 170	0.110 11.0
IgG ABX Pentra IgG CP (Ref. A11A01924)	•			g/L mg/dL	6.98 698	5.58 – 8.38 558 – 838	0.467 46.7
		•		g/L mg/dL	6.98 698	5.58 – 8.38 558 – 838	0.467 46.7
			•	g/L mg/dL	6.98 698	5.58 – 8.38 558 – 838	0.467 46.7
				g/L mg/dL	6.98 698	5.58 – 8.38 558 – 838	0.467 46.7
IgM ABX Pentra IgM CP (Ref. A11A01925)	•			g/L mg/dL	0.80 80	0.58 – 1.02 58 - 102	0.073 7.3
		•		g/L mg/dL	0.80 80	0.58 – 1.02 58 - 102	0.073 7.3
			•	g/L mg/dL	0.80 80	0.58 – 1.02 58 - 102	0.073 7.3
				g/L mg/dL	0.80 80	0.58 – 1.02 58 - 102	0.073 7.3
Transferrin ABX Pentra Transferrin CP (Ref. A11A01926)	•			g/L mg/dL	1.27 127	1.02 – 1.52 102 - 152	0.083 8.3
		•		g/L mg/dL	1.27 127	1.02 – 1.52 102 - 152	0.083 8.3
			•	g/L mg/dL	1.27 127	1.02 – 1.52 102 - 152	0.083 8.3
				g/L mg/dL	1.27 127	1.02 – 1.52 102 - 152	0.083 8.3