

ABX Pentra Glucose PAP CP	A91A00588BEN	11/10/2006
----------------------------------	---------------------	-------------------

1. Product and company identification	
1.1. Identification of the product	
Product name:	ABX Pentra Glucose PAP CP
Product code:	Ref.ABX: A11A01668
1.2. Company identification	
	HORIBA ABX - Rue du Caducée - Parc Euromédecine 34184 MONTPELLIER CEDEX 4 - FRANCE Tel: (33) 4 67 14 15 16 Fax: (33) 4 67 14 15 17
1.3. Emergency phone number	
	Contact the nearest first-aid station.

2. Composition/information on ingredients	
2.1. Description of the kit	
	REAGENT Mono-reagent cassette; 1 x 99 ml; Phosphate buffer (13.8 mmol/l), Phenol (10 mmol/l), 4-amino-antipyrine (0.3 mmol/l), Glucose oxydase ($\geq 10\ 000$ U/l), Peroxydase (≥ 700 U/l)
2.2. Hazardous ingredients (name, concentration)	
Sodium azide	Less than 0.1%
2.3. Dangerous preparation according to 67/548/EEC - 99/45/EEC (Yes/No): No	

3. Hazards identification	
Most important hazards:	This reagent is particularly stable and no toxic in normal conditions of use.
Specific risks:	<u>For man, effect of acute exposure:</u> Eyes, skin contact: May cause irritation. Inhalation: May cause irritation, may be harmful. Ingestion: May be harmful. <u>For environment:</u> None hazard determined.

4. First-aid-measures	
Inhalation:	Provide fresh air. If breathing is difficult, get medical attention.
Skin contact:	Immediately wash skin with plenty of water.
Eye contact:	Immediately wash eyes, also under the eyelids, with water for at least 15 minutes. Get medical attention.
Ingestion:	Get medical attention immediately.

ABX Pentra Glucose PAP CP	A91A00588BEN	11/10/2006
----------------------------------	---------------------	-------------------

5. Fire-fighting measures

Extinguishing media:	Water spray or foam.
Specific risks:	Not determined.
Special protective equipment:	Not determined.
Additional recommendations:	Thermal decomposition: Dangerous decomposition products and gases which are damaging to health, in large quantities, are not anticipated.

6. Accidental release measures

Personal precautions:	Avoid eyes, skin or clothes contact. Rinse well with water. Provide fresh air.
Environmental precautions:	Not determined.
Methods for cleaning/absorption:	Actions to be taken on spillage: soak up and remove with absorbent material. Absorbent material: no special requirements. Disposal measure: collect waste in closed disposal containers and dispose of according to local regulations.

7. Handling and storage

7.1. Handling

Technical measures:	Not determined.
Precautions:	Not determined.
Handling recommendations:	Avoid eyes, skin contact. Wash hands before and after use.

7.2. Storage

Technical measures:	Not determined.
Storage conditions:	Store protected from light and tightly closed at 2-8°C.
Incompatible materials:	Not determined.
Packaging materials:	Not determined.

8. Exposure controls/personal protection

Specific control parameters:	Not determined.
Technical measures:	Not determined.
Respiratory protection:	Have good local ventilation, wear mask or appropriate respirator.
Hand protection:	Protective gloves.
Eye protection:	Safety goggles.
Skin protection:	Lab coat.
Hygiene measures:	Wash hands before and after use.

ABX Pentra Glucose PAP CP	A91A00588BEN	11/10/2006
----------------------------------	---------------------	-------------------

9. Physical and chemical properties			
Physical state:	Liquid	Autoflammability:	Not determined
Colour:	Colourless	Explosive properties:	Not determined
Odour:	Not determined	Oxidizing properties:	Not determined
pH:	Not determined	Vapour pressure:	Not determined
Boiling point:	Not determined	Relative density:	Not determined
Melting point:	Not determined	Solubility:	Not determined
Flashing point:	Not determined	Partition coefficient:	Not determined
Flammability (solid, gaz):	Not determined	Other data:	None relevant for the safety.

10. Stability and reactivity	
Conditions to avoid:	The product is stable if stored at 2-8°C, avoid exposure to heat.
Materials to avoid:	Not determined.
Hazardous decomposition products:	Not determined.
Other data:	Stability: The product is stable until the expiry date on the label. No hazardous reactions under normal use.

11. Toxicological information	
Acute toxicity:	Not determined.
Chronic toxicity:	Not determined.
Inhalation:	The reagent contains small amount of phenol (< 0.1%). May be harmful. The reagent contains small amount of sodium azide (< 0.1%) as a preservative and may be harmful if swallowed.
Skin contact:	The reagent contains small amount of phenol (< 0.1%). Contact of large amount may be harmful. Handling may cause irritation of skin and mucous membranes.
Eye contact:	Handling may cause irritation of eyes.
Ingestion:	The reagent contains small amount of phenol (< 0.1%). May be harmful. The reagent contains small amount of sodium azide (< 0.1%) as a preservative. Ingestion of large amount may cause nausea, vomiting and in certain circumstances respiratory difficulties, high pulse rate and/or hypersensitivity.
Other data:	Toxicological data for this reagent are not available.

12. Ecological information	
Mobility:	Not determined.
Degradability:	Not determined.
Bioaccumulative potential:	Not determined.
Ecotoxicology:	When handled and used properly, no ecological problems are to be expected.
Other hazardous effects:	Ecological data for this product are not available. Disposal procedure have to be respected (see §13).

ABX Pentra Glucose PAP CP	A91A00588BEN	11/10/2006
----------------------------------	---------------------	-------------------

13. Disposal considerations

Surplus or waste (residues):	This reagent contains less than 0.1% sodium azide as a preservative. As sodium azide may react with lead and copper to form explosive metal azides, this reagent should be disposed of by flushing with copious amounts of water. The product must be disposed of as a laboratory chemical according to local regulations.
Contaminated packaging:	Type of primary container: The reagent is packaged in plastic cassette (composed of polyethylene high density).

14. Transport information (International regulations)

General information:	This product is not subject to current regulations for transportation of hazardous goods.
By air (IATA):	Not applicable.
By land (ADR):	Not applicable.
By sea (IMDG):	Not applicable.

15. Regulatory information

Symbols:	None	
R phrases:	None	
S phrases:	None	
Substances:	None	
Other regulatory requirements:	In no way does this information exempt the user from knowing about or applying all the national or international regulations relating to his/her activity.	

16. Other information

This sheet is a complement to instructions for use but does not replace them. Information contained herewith is based on our current knowledge of the product, at the date shown, and is correct to the best of our knowledge. Furthermore, the user's attention is drawn to the dangers of using this product for anything other than its intended use. The user must accept the sole responsibility and take precautions accordingly for the use of this product.