

| | | |
|---------------------------|---------------------|--------------------------------|
| ABX MINOTROL RETIC | A91A00286DEN | Revision 08/09/2008 |
|---------------------------|---------------------|--------------------------------|

| | |
|--|---|
| 1. Product and company identification | |
| 1.1. Identification of the product | |
| Product name: | ABX MINOTROL RETIC |
| Product code: | 2072201: Level 2, Twin-pack (2 x 3ml vial + value sheet) 2072202: Level 1 + Level 3, Twin-pack (2 x 3ml vial + value sheet) 2072001: Level 1 (1 x 3ml vial) 2072002: Level 2 (1 x 3ml vial) 2072003: Level 3 (1 x 3ml vial) |
| 1.2. Use of the product | |
| ABX Minotrol retic is a tri-level, assayed hematology control designed to document and monitor values obtained from manual and automated | |
| 1.3. Company identification | |
| | HORIBA ABX - Rue du Caducée - Parc Euromédecine 34184 MONTPELLIER CEDEX 4 - FRANCE Tel: +33 (0) 4 67 14 15 16 Fax: +33 (0) 4 67 14 15 17 Email: documentation@fr.abx.fr |
| 1.4. Emergency phone number | |
| | Contact the nearest first-aid station |

| | |
|---|--|
| 2. Hazards identification | |
| Dangerous preparation according to 67/548/EEC - 99/45/EEC (Yes/No): No | |
| Most important hazards: | Non applicable. |
| Specific risks | Potentially biohazardous material. Handle as if potentially infectious. |

| 3. Composition/information on ingredients | | | | | |
|---|---|---------------|--------|------|---------------------|
| 3.1. Description of the kit | | | | | |
| | Contains human erythrocytes and mammalian erythrocytes suspended in a plasma-like fluid with preservatives. «1» : Low level, «2» : Medium level, «3» : High level. | | | | |
| 3.2. Hazardous ingredients (name, concentration) | | | | | |
| CAS N. | Name | Concentration | Symbol | Risk | Registration number |
| To the best of our knowledge, this product contains no hazardous ingredients as defined by the OSHA Hazard Communication Standard, 29 CFR 1910.1200 and by EU Directive 67/548/EEC. | | | | | |

| | | |
|---------------------------|---------------------|--------------------------------|
| ABX MINOTROL RETIC | A91A00286DEN | Revision 08/09/2008 |
|---------------------------|---------------------|--------------------------------|

4. First-aid-measures

| | |
|---------------|--|
| Inhalation: | Not determined. |
| Skin contact: | Immediately wash skin with soap and water. |
| Eye contact: | Immediately wash eyes, also under the eyelids, with water for at least 15 minutes. |
| Ingestion: | Rinse mouth out with water if the victim is conscious, otherwise call doctor. |

5. Fire-fighting measures :

| | |
|-------------------------------|--------------------------------|
| Extinguishing media: | This product is non-flammable. |
| Specific risks: | None determined. |
| Special protective equipment: | None determined. |
| Additional recommendations: | None. |

6. Accidental release measures

| | |
|----------------------------------|---|
| Personal precautions: | Wear laboratory gloves and protective goggles. |
| Environmental precautions: | Neutralise product before discarding. |
| Methods for cleaning/absorption: | After neutralisation, absorb the product and rinse the whole area with water. |

7. Handling and storage

7.1. Handling

| | |
|---------------------------|---|
| Technical measures: | Avoid contact with the skin, eyes. Laboratory gloves and protective goggles. |
| Precautions: | Always handle the product in a well-ventilated area. |
| Handling recommendations: | Not determined. |

7.2. Storage

| | |
|-------------------------|--|
| Technical measures: | Not determined. |
| Storage conditions: | To maintain the quality of the product, store ABX MINOTROL RETIC upright at 2° to 8°C when not in use. Protect vials from overheating and freezing. Unopened vials are stable until the expiration date. Opened tubes are stable for 16 sampling events over a maximum of 16 days ^a . |
| Incompatible materials: | Not determined. |
| Packaging materials: | Not determined. |

a. Modification from index C to D : stability changed.

| | | |
|---------------------------|---------------------|--------------------------------|
| ABX MINOTROL RETIC | A91A00286DEN | Revision 08/09/2008 |
|---------------------------|---------------------|--------------------------------|

| | |
|---|---|
| 8. Exposure controls/personal protection | |
| 8.1. Exposure limit values | |
| Specific control parameters: | Not determined. |
| 8.2. Exposure controls | |
| 8.2.1. Occupational exposure controls | |
| Technical measures: | ABX MINOTROL RETIC should be handled in the same manner as a human blood sample. Eye wash, safety shower Laboratory gloves, goggles in case of prolonged contact. |
| Respiratory protection: | Not determined. |
| Hand protection: | Protective gloves. |
| Eye protection: | Safety goggles. |
| Skin protection: | Lab coat. |
| Hygiene measures: | Follow good laboratory practices. |
| 8.2.2. Environmental exposure controls | |
| Not determined. | |

| | | | |
|--|-----------------|------------------------|----------------------|
| 9. Physical and chemical properties | | | |
| 9.1. General information | | | |
| Physical state: | Liquid. | | |
| Colour: | Red/Black | Odour: | Odourless |
| 9.2. Important health, safety and environmental information | | | |
| pH: | Not determined. | Relative density: | Not determined. |
| Boiling point: | Not determined. | Solubility: | Miscible with water. |
| Flash point: | Not determined. | Partition coefficient: | Not determined. |
| Flammability (solid, gaz): | Non-flammable. | Viscosity: | Not determined. |
| Explosive properties: | Not determined. | Vapour density: | Not determined. |
| Oxidizing properties: | Not determined. | Evaporation rate: | Not determined. |
| 9.3. Other information | | | |
| Not determined. | | | |

| | | |
|---------------------------|---------------------|--------------------------------|
| ABX MINOTROL RETIC | A91A00286DEN | Revision 08/09/2008 |
|---------------------------|---------------------|--------------------------------|

| | |
|-------------------------------------|-----------------|
| 10. Stability and reactivity | |
| Conditions to avoid: | Expiry date. |
| Materials to avoid: | Not determined. |
| Hazardous decomposition products: | Not classified. |
| Other data: | None. |

| | |
|--------------------------------------|--|
| 11. Toxicological information | |
| Acute toxicity: | Not classified according to directive 99/45/EEC. |
| Chronic toxicity: | Not classified according to directive 99/45/EEC. |
| Inhalation: | Not determined. |
| Skin contact: | Not determined. |
| Eye contact: | Not determined. |
| Ingestion: | Not determined. |
| Other data: | None. |

| | |
|---|----------------------|
| 12. Ecological information | |
| 12.1. Ecotoxicity: | Not determined. |
| 12.2. Mobility: | Not determined. |
| 12.3. Degradability: | Biodegradable > 90%. |
| 12.4. Bioaccumulative potential: | Not determined. |
| 12.5. Results of PBT assessment: | None. |
| 12.6. Other hazardous effects: | Not determined. |

| | |
|------------------------------------|--|
| 13. Disposal considerations | |
| Surplus or waste (residues): | The product must be disposed of as a laboratory chemical according to local regulations. |
| Contaminated packaging: | Incineration. |

| | | |
|---------------------------|---------------------|--------------------------------|
| ABX MINOTROL RETIC | A91A00286DEN | Revision 08/09/2008 |
|---------------------------|---------------------|--------------------------------|

| | |
|--|-----------------|
| 14. Transport information (International regulations) | |
| General information: | None. |
| By air (IATA): | Not applicable. |
| By land (ADR): | Not applicable. |
| By sea (IMDG): | Not applicable. |

| | |
|---------------------------------------|--|
| 15. Regulatory information | |
| Symbols: | None. |
| R phrases: | None. |
| S phrases: | None. |
| Substances: | None. |
| Specific provisions and related laws: | None. |
| Other regulatory requirements: | In no way does this information exempt the user from knowing about or applying all the national or international regulations relating to his/her activity. |

| | |
|---|--|
| 16. Other information | |
| 16.1. List of relevant R phrases | |
| None. | |
| 16.2. Recommended restriction on use | |
| See specific information given in the product user information. | |
| 16.3. Further information | |
| <p>This sheet is a complement to instructions for use but does not replace them. Information contained herewith is based on our current knowledge of the product, at the date shown, and is correct to the best of our knowledge.</p> <p>Furthermore, the user's attention is drawn to the dangers of using this product for anything other than its intended use. The user must accept the sole responsibility and take precautions accordingly for the use of this product.</p> | |
| 16.4. Sources of key data | |
| This document has been compiled with information extracted from documents given by our raw material suppliers. | |
| 16.5. Revision | |
| The sections that have been modified from the previous version are indicated by a vertical line on the left side of the table. | |