

ABX Miniclean 1L	A91A00204BEN	16/12/2005
-------------------------	---------------------	-------------------

1. Product and company identification	
1.1. Identification of the product	
Product name:	ABX Miniclean 1L
Product code:	Ref.ABX: 0403010
1.2. Company identification	
	HORIBA ABX - Rue du Caducée - Parc Euromédecine 34184 MONTPELLIER CEDEX 4 - FRANCE Tel: (33) 4 67 14 15 16 Fax: (33) 4 67 14 15 17
1.3. Emergency phone number	
	Contact the nearest first-aid station

2. Composition/information on ingredients	
2.1. Description of the kit	
	[REAGENT]: 1 x 1L Cleaning solution : Organic buffer <20%, Proteolytic enzyme <1%
2.2. Hazardous ingredients (name, concentration)	
2.3. Dangerous preparation according to 67/548/EEC - 99/45/EEC (Yes/No): No	

3. Hazards identification	
Most important hazards:	None.
Specific risks	Product may be harmful if swallowed and can induce vomiting, fainting and diarrhoea.

4. First-aid-measures	
Inhalation:	No specific requirements.
Skin contact:	Wash skin with water and soap.
Eye contact:	Immediately wash eyes, also under the eyelids, with water for at least 15 minutes.
Ingestion:	If the victim is conscious, rinse out mouth with water. Get medical attention in case of fainting.

ABX Miniclean 1L	A91A00204BEN	16/12/2005
-------------------------	---------------------	-------------------

5. Fire-fighting measures : Product non flammable

Extinguishing media:	Use whatever is required for the surrounding area.
Specific risks:	None determined.
Special protective equipment:	None.
Additional recommendations:	None.

6. Accidental release measures

Personal precautions:	Wear a lab coat, protective gloves and safety goggles.
Environmental precautions:	Not determined.
Methods for cleaning/absorption:	Soak up the product and put in the waste. Wash the affected area with water.

7. Handling and storage

7.1. Handling

Technical measures:	None determined.
Precautions:	The usual laboratory precautions should be observed.
Handling recommendations:	Wear gloves and protective goggles in case of prolonged handling.

7.2. Storage

Technical measures:	No special requirements.
Storage conditions:	Store at 18-25°C.
Incompatible materials:	Not determined.
Packaging materials:	Not determined.

8. Exposure controls/personal protection

Specific control parameters:	None determined.
Technical measures:	Use general room ventilation.
Respiratory protection:	None required.
Hand protection:	Protective gloves.
Eye protection:	Safety goggles.
Skin protection:	Lab coat.
Hygiene measures:	Wash hands before and after use.

ABX Miniclean 1L	A91A00204BEN	16/12/2005
-------------------------	---------------------	-------------------

9. Physical and chemical properties

Physical state:	Liquid	Autoflammability:	Not determined
Colour:	Colourless	Explosive properties:	Not determined
Odour:	Smells of detergent	Oxidizing properties:	Not determined
pH:	9.6 (+/- 0.4)	Vapour pressure:	Not determined
Boiling point:	Not determined	Relative density:	Not determined
Melting point:	Not determined	Solubility:	Aqueous solution miscible with water
Flashing point:	Not determined	Partition coefficient:	Not determined
Flammability (solid, gaz):	Non-flammable	Other data:	Not determined

10. Stability and reactivity

Conditions to avoid:	None determined.
Materials to avoid:	None determined.
Hazardous decomposition products:	None determined.
Other data:	None.

11. Toxicological information

Acute toxicity:	Not determined.
Chronic toxicity:	Not determined.
Inhalation:	Not determined.
Skin contact:	Not determined.
Eye contact:	Not determined.
Ingestion:	Not determined.
Other data:	None.

12. Ecological information

Mobility:	Not determined.
Degradability:	Biodegradable >90%
Bioaccumulative potential:	Not determined.
Ecotoxicology:	Not determined.
Other hazardous effects:	None determined.

ABX Miniclean 1L	A91A00204BEN	16/12/2005
-------------------------	---------------------	-------------------

13. Disposal considerations	
Surplus or waste (residues):	The product must be disposed of as a laboratory chemical according to local regulations.
Contaminated packaging:	Incineration.

14. Transport information (International regulations)	
General information:	None.
By air (IATA):	Not applicable.
By land (ADR):	Not applicable.
By sea (IMDG):	Not applicable.

15. Regulatory information		
Symbols:	None	
R phrases:	None	
S phrases:	None	
Substances:		
Other regulatory requirements:	In no way does this information exempt the user from knowing about or applying all the national or international regulations relating to his/her activity.	

16. Other information	
<p>This sheet is a complement to instructions for use but does not replace them. Information contained herewith is base on our current knowledge of the product, at the date shown, and is correct to the best of our knowledge.</p> <p>Furthermore, the user's attention is drawn to the dangers of using this product for anything other than its intended use.</p> <p>The user must accept the sole responsibility and take precautions accordingly for the use of this product.</p>	