



HORIBA ABX SAS
 Parc Euromédecine - Rue du Caducée
 B.P. 7290
 34184 Montpellier Cedex 4
 FRANCE
 www.horiba.com



LOT L: 22564

REF : A11A01928

04710819D

SAP 1300023129



2019/05/31

(FR) Constituant (EN) Constituent (DE) Bestandteil (ES) Componente (IT) Componente (PT) Componente (DA) Bestanddel (SV) Konstituent (EL) συστατικό (PL) Składnik (NO) Bestanddel	ABX Pentra 400	Pentra C400	Pentra C200	Unité Unit Einheit Unidad Unità Unidade Enhed Enhet Μονάδα Jednostka Enhet	Valeur titrée Assigned value Deklariertes Wert Valor asignado Valori assegnati Valor determinado Tildelt værdi Tildelat värde Αναγραφόμενη τιμή Wartość przypisana Tildelt verd	Intervalle de confiance Confidence range Vertrauensintervalle Intervalos de confianza Livelli di fiducia Intervalo de confiança Konfidensområde Konfidensintervall Εύρος εμπιστοσύνης Zakres ufności Konfidensspekter	
C				U	\oplus	→ ±3 SD ←	1 SD
IgA ABX Pentra IgA CP (Ref. A11A01923)	•			g/L 1.31	1.31	0.99 – 1.62	0.107
		•		mg/dL 131	131	99 - 162	10.7
			•	g/L 1.31	1.31	0.99 – 1.62	0.107
			•	mg/dL 131	131	99 - 162	10.7
IgG ABX Pentra IgG CP (Ref. A11A01924)	•			g/L 7.35	7.35	5.88 – 8.81	0.490
		•		mg/dL 735	735	588 - 881	49.0
			•	g/L 7.35	7.35	5.88 – 8.81	0.490
			•	mg/dL 735	735	588 - 881	49.0
IgM ABX Pentra IgM CP (Ref. A11A01925)	•			g/L 1.01	1.01	0.73 – 1.29	0.093
		•		mg/dL 101	101	73 - 129	9.3
			•	g/L 1.01	1.01	0.73 – 1.29	0.093
			•	mg/dL 101	101	73 - 129	9.3
Transferrin ABX Pentra Transferrin CP (Ref. A11A01926)	•			g/L 1.33	1.33	1.06 – 1.60	0.090
		•		mg/dL 133	133	106 - 160	9.0
			•	g/L 1.33	1.33	1.06 – 1.60	0.090
			•	mg/dL 133	133	106 - 160	9.0