



**HORIBA ABX SAS**  
 Parc Euromédecine - Rue du  
 Caducée  
 B.P. 7290  
 34184 Montpellier Cedex 4  
 FRANCE  
 www.horiba.com



**LOT** 1: 1700801 (rev.01)



1 7 0 0 8 0 1

**REF** :

1300026187



2018/31/10

1300026587



1300026587

(FR) Constituant (EN) Constituent (DE) Bestandteil (ES) Componente (IT) Componente (PT) Componente (DA) Bestanddel (SV) Konstituent (EL) συστατικό (PL) Składnik (NO) Bestanddel	ABX Pentra 400	Pentra C400	Pentra C200	Unité	Valeur titrée	Intervalle de confiance	
				Unit	Assigned value	Confidence range	1 SD
				Einheit	Deklariertes Wert	Vertrauensintervalle	
				Unidad	Valor asignado	Intervalos de confianza	
				Unità	Valori assegnati	Livelli di fiducia	
				Unidade	Valor determinado	Intervalo de confiança	
				Enhed	Tildelt værdi	Konfidensområde	
				Enhet	Tildelat värde	Konfidensintervall	
				Μονάδα	Αναγραφόμενη τιμή	Εύρος εμπιστοσύνης	
				Jednostka	Wartość przypisana	Zakres ufności	
				Enhet	Tildelt verd	Konfidensspekter	
C				U	$\oplus$	Confidence range	1 SD
CK (IFCC/DGKC) (37°C)	•			U/L µkat/L	<b>153</b> <b>2.56</b>	138 - 168 2.32 - 2.80	5 0.08
		•		U/L µkat/L	<b>153</b> <b>2.56</b>	138 - 168 2.32 - 2.80	5 0.08
			•	U/L µkat/L	<b>155</b> <b>2.59</b>	140 - 170 2.35 - 2.83	5 0.08
CK-MB (Immunology method) (37°C)	•			U/L µkat/L	<b>32.9</b> <b>0.55</b>	27.5 - 38.3 0.46 - 0.64	1.8 0.03
		•		U/L µkat/L	<b>32.9</b> <b>0.55</b>	27.5 - 38.3 0.46 - 0.64	1.8 0.03
			•	U/L µkat/L	<b>40.1</b> <b>0.67</b>	32.9 - 47.3 0.55 - 0.79	2.4 0.04